



April 21, 2023

Dräger Medical Verification of LANCOM LX-6402 AP with LANCOM WLC-1000 Controller for OneNet® (shared Dräger Infinity network)

Dear Madam/Sir,

Dräger Medical Systems is pleased to announce the verification of LANCOM LX-6402 AP with LANCOM WLC-1000 controller running the below mentioned firmware versions. These components have been extended to the Infinity OneNet approved network equipment checklist, OneNet is the award-winning solution introduced by Dräger Medical that enables the integration of Infinity patient monitors as part of a single hospital wide network. The following network components are verified for use with the Dräger Medical Infinity Network:

WLAN Controller	Access Point	Firmware Version	Comments
LANCOM Controller WLC-1000	LANCOM LX-6402 AP	Controller SW version : 10.72.0015 AP SW Version : 6.10.0030	Passed.

The goal of this verification initiative is to investigate situations where customers want to operate Dräger and non-Dräger clients on shared IEEE 802.11 a/b/g/n/ac/ax access points. The above referenced access point has passed the simulated tests where an access points were loaded with eight wireless Infinity Delta patient monitors, 16 Infinity M300 monitors, eight M540 standalone monitors, as well as competing traffic. Please be advised it remains the customer's responsibility to actively monitor potential sources of wireless interference as well as changes to wireless traffic loads.

Yours sincerely,

Sabari

Srinivasan

Digitally signed by
Sabari Srinivasan
Date: 2023.04.21
22:35:19 -04'00'

Sabari Gopalakrishnan Srinivasan
Sr Product Manager
Business Unit Patient Monitoring
Mobile : 978 409 0643
Email : sabari.srinivasan@draeger.com

Dräger Medical Systems, Inc.

6 Tech Drive
Andover, MA 01810 USA
Phone: (978) 379-8000
www.draeger.com